

# healthcare

# HTM64 (SU H) firth scrub-up trough

SU H scrub-up trough in units of 80cm, 160cm and 240cm for 1, 2 or 3 persons, complete with 1½in strainer waste and trap cover, with earthing terminal.

Scrub-up trough assembly

Outlet handed left or right

Screw to wall

Stainless steel



# Illustrated **S2860 (SU H 3)** 240cm

with right hand outlet, strainer waste & trap cover

## A4555AA (TB H6)

Markwik 21, integral thermostat, bib mixer, time flow sensor, single flow horizontal outlet, 2 x 1/2 in inlets

**S8925 (TRR2/P)** 1½in plastic resealing bottle trap with removable sump, 75mm seal

### Weight

Stainless steel type 1.4301 (304) Trough 1.5mm

# Weight

240cm trough 46.0kg 160cm trough 34.0kg 80cm trough 16.0kg

Exposed surfaces polished with 240 grain

### **Standards**

Stainless steel conforms to BS EN 10088:Part 2:2005

### **Options**

**S2858 (SU H 3)** 240cm with left hand outlet, strainer waste & trap cover **S2854 (SU H 2)** 160cm with left hand outlet, strainer waste & trap cover **S2856 (SU H 2)** 160cm with right hand outlet. strainer waste & trap cover **S2850 (SU H 1)** 80cm with left hand outlet, strainer waste & trap cover S2852 (SU H 1) 80cm with right hand outlet, strainer waste & trap cover

# **Special notes**

For use in connection with surgical washing of forearms and hands see page 19:1:16C for TB H6 details









